



# SERIES 5 DENTAL CHAIR

## Technical Specification

**DCI** edge



CHAIR WEIGHT: 300lbs / 136.4kg  
UPHOLSTERY WEIGHT: 30lbs / 13.6kg  
SHIPPING WEIGHT:

POWER REQUIREMENTS: 115VAC, 8A, 50-60Hz  
230VAC, 2.5A, 50-60Hz  
POWER CORD LENGTH: 5 Feet (1.5m)

IEC MEDICAL CLASSIFICATION I  
TYPE: B Chair upholstery  
OPERATION: Intermittent  
DUTY CYCLE: 20 Seconds ON, 5 Minutes OFF  
SPLASH PROTECTION: IPX0 (FOOT SWITCH IPX1)

CERTIFICATIONS: IEC 60601-1 3rd Edition,  
CAN/CSA C22.2 No. 601.1

