

DCI EDGE



DCI EDGE SERIES 4 CABINETRY

CABINETS & DELIVERY

12 O'CLOCK CABINET

The Series 4 offers affordable quality with dental specific cabinetry at a contractor's price. DCI's 12 o'clock cabinets feature simplicity, functionality and the professional aesthetic our customers have come to expect when ordering the DCI Edge line.

12 O'CLOCK DELIVERY

The DCI Edge 12 o'clock delivery system provides a sturdy platform for the dental team to operate.

Heavy Duty Arm w/ Worksurface

High quality construction with reliable pneumatic components will stand the test of time.

Integrated Touchpad

Provides ergonomic access to easily control your equipment, eliminating unnecessary movement.

Assistant's Work Surface & Instrumentation

Provides a sturdy and efficient work surface for the Dental Assistant. Integrated quick disconnect and electrical outlets make adding mobile accessories a breeze.



SIDE CABINETS

DCI Edge Series 4 dental cabinets were designed to optimize function and flow of the dental operator. Unlike cabinets built by contractors, DCI Edge cabinets are serialized and carry the same tax benefits as your other operatory equipment.

Constructed with a structural steel sub-base, high pressure industrial grade laminate and 1mm PVC edgebanding, these cabinets will look great for years to come, all while maintaining an affordable price.





Started 35 years ago, DCI was founded on three guiding principles: making our customers' lives easier, providing dental care for those in need, and enhancing the lives of our employees. These principles give us our edge. It's an edge we want to pass along to you.


DCI is a family owned company for two generations which is based in Newberg, Oregon, the heart of dental equipment manufacturing in America. Over the years, John & Jason Spencer redefined the dental landscape by designing dental equipment that challenged the status quo. Today, the family brings their design expertise and commitment to value and service to the forefront with the DCI Edge Cabinetry.



Jason, Amy and John Spencer

Cabinet Color Selections

Standard Wilsonart® Laminate Selections for Vertical Surfaces



Standard Wilsonart® Solid Surface Selections for Countertops

